

## Technical description

# INDIVIDUAL HIP®

### Description

Individualized femoral stem, delivered together with an individualized smooth compactor or rasp.



### References

#### 3000 0201 INDIVIDUAL HIP® TiHA

Individualized cementless stem  
Titanium alloy (Ti6Al4V-ISO 5832-3), coated with porous titanium and hydroxyapatite  
12/14 5°40' taper

#### Design options

##### Coating

- Partial HA coating
- Full HA coating

##### Collar

- No collar
- Normal collar



#### 3000 0301 INDIVIDUAL HIP® Cemented

Individualized cemented stem  
Stainless steel (M30NW-ISO 5832-9)  
12/14 5°40' taper

#### Design options

##### Collar

- No collar
- Normal collar



## Instrumentation

According to the surgeon preferences in terms of instrumentation and surgical approach, the INDIVIDUAL HIP® stem is delivered with one of the following model of compactor or rasp.

**3000 0501 INDIVIDUAL HIP® Compactor**

Individualized smooth compactor  
Stainless steel (316L-NF S94-090)  
12/14 5°40' taper



**3000 0701 INDIVIDUAL HIP® Rasp**

Individualized piqued rasp  
Stainless steel (316L-NF S94-090)  
12/14 5°40' taper



**3000 0511 INDIVIDUAL HIP® Compactor NT**

Individualized smooth compactor, without trial neck  
Stainless steel (316L-NF S94-090)



**3000 0711 INDIVIDUAL HIP® Rasp NT**

Individualized piqued rasp, without trial neck  
Stainless steel (316L-NF S94-090)



Moreover, the following instrumentation is required for the implantation of an INDIVIDUAL HIP® stem:

- **7011 0000** INDIVIDUAL HIP® instrumentation



7011 0000 INDIVIDUAL HIP® instrumentation