


Symbios Connect Tutorial

Base functions



First login

How to proceed to the first login

Enter your login with blank password

1

Type here the e-mail address registered for your Symbios Connect account



English

Symbios Connect

Connection

john.doe@symbios.ch

Password

Connection

 Remain connected[Forgot password](#)

2

Leave the password as blank and click on the 'Connection' button

Check for the verification code on your mailbox

3

An e-mail has been sent on your e-mail account. **Copy the verification code.**

SymbiosConnect - Verification code



Symbios Connect <connect@symbios.ch>

À: sc validation1

Aujourd'hui à 14:05

Hello,

The verification code required to setup or reset your password is: **078JqT.**

If you are not modifying your SymbiosConnect account, please contact our technical team to ensure your account is safe.

Sincerely,

Symbios Connect Customer Service



Set your password and enter the verification code

4 Fill the form with your password and confirm the password

The screenshot shows the Symbios Connect login page. The header includes the Symbios logo and the text 'Symbios Connect'. A callout box labeled '4' points to the password confirmation form. The form has the following fields and elements:


- Password confirmation** header
- Requirement text: "Your password should contain at least 8 characters, one number, and UPPER and lower case letters"
- Password** input field with a strength indicator: "8 characters Numbers UPPER / lower case"
- Identical password** input field
- Code sent by email** input field
- Confirm** button

5 Paste the verification code in the textbox and click on the 'Confirm' button

Finalize by logging in with your new credentials

6

Now you can log on the login page again



The screenshot shows the Symbios Connect login page. In the top left corner, there is the Symbios logo with the tagline "custom-made for you". The main heading "Symbios Connect" is prominently displayed on the left side. On the right side, there is a message indicating a successful password change: "Change password" followed by "Password changed successfully." and a green link that says "Return to the sign-in page." A callout box with the number 6 points to the top of the page, containing the text "Now you can log on the login page again".

User interface overview

First steps to understand the interface

List of patient cases

Symbios Connect Dr CHARNLEY John ▾

Surgical cases

 **List of patients for which a surgery is planned and/or a custom product ordered**

 New case 

Date quick filter

« July 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
27	28	29	30	1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31
1	2	3	4	5	6	7

Today

Status quick filters

Active cases (18)

- Pending approvals (4)
- In process (13)
- Delayed (1)

Inactive cases (1)

- Operated (1)
- Cancelled (0)
- Suspended (0)
- To complete (0)

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN®	Shipped
	HARI Mata	KNEE	R	ORIGIN®	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN®	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN®	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN®	Planning approved
	CHANEL Coco	HIP	L	INDIVIDUAL HIP™	Planning in progress
	MONROE Marilyn	KNEE	R	ORIGIN®	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN®	Pending planning approval
	MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available
15.09.2022	CHURCHILL Winston	HIP	R	INDIVIDUAL HIP™	Pending planning approval
	FITZGERALD Ela	HIP	R	HIP-PLAN	Waiting for CT scan
04.10.2022	KELLY Grace	KNEE	L	ORIGIN®	Pending planning approval

Rows / page: **15** 30 50

< Previous 1 2 Next >

Filter and search functions

Symbios Connect Dr CHARNLEY John ▾

filters **Search textbox**

Surgical cases New case

Date quick filter

« July 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
27	28	29	30	1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31
1	2	3	4	5	6	7

Today

Status quick filters

Active cases [18]

- Pending approvals (4)
- In process (13)
- Delayed (1)

Inactive cases [1]

- Operated (1)
- Cancelled (0)
- Suspended (0)
- To complete (0)

Rows / page: 15 30 50

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN®	Shipped
	HARI Mata	KNEE	R	ORIGIN®	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN®	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN®	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN®	Planning approved
	CHANEL Coco	HIP	L	INDIVIDUAL HIP™	Planning in progress
	MONROE Marilyn	KNEE	R	ORIGIN®	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN®	Pending planning approval
	MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available
15.09.2022	CHURCHILL Winston	HIP	R	INDIVIDUAL HIP™	Pending planning approval
	FITZGERALD Ela	HIP	R	HIP-PLAN	Waiting for CT scan
04.10.2022	KELLY Grace	KNEE	L	ORIGIN®	Pending planning approval

< Previous 1 2 Next >

User profile menus / settings

Symbios Connect

Surgical cases

 **Click here to display user profile menus + settings**

Pr CHARNLEY John ▾

- Pr CHARNLEY John
marketing@symbios.ch
- Preferences
- Downloads
- Help
- Terms of Use
- Sign out

Search

 New case

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]
- Suspended [0]

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN®	Shipped
	HARI Mata	KNEE	R	ORIGIN®	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN®	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN®	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN®	Planning approved
	MONROE Marilyn	KNEE	R	ORIGIN®	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN®	Pending planning approval
	MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available

< Previous 1 2 Next >

Quick view on an order

Symbios Connect Dr CHARNLEY John ▾

Place mouse over the patient's row to display case quick view

Surgical cases

New case

Date quick filter

« July 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
27	28	29	30	1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31
1	2	3	4	5	6	7

Today

Status quick filters

Active cases (18)

- Pending approvals (4)
- In process (13)
- Delayed (1)

Inactive cases (1)

- Operated (1)
- Cancelled (0)
- Suspended (0)
- To complete (0)

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri				Shipped
	HARI Mata				Waiting for CT scan
30.08.2022	DESCARTES René				Manufacturing in progress
	BAKER Joséphine				Ready for delivery
	CHOPIN Frédéric				CT scan available
	ROOSEVELT Franklin				CT scan available
06.09.2022	JEFFERSON Thomas				Planning approved
	CHANEL Coco				Planning in progress
	MONROE Marilyn				Pending planning approval
13.09.2022	EINSTEIN Albert				Pending planning approval
	MANDELA Nelson				CT scan available
15.09.2022	CHURCHILL Winston				Pending planning approval
	FITZGERALD Eta	HIP	R	HIP-PLAN	Waiting for CT scan
04.10.2022	KELLY Grace	KNEE	L	ORIGIN®	Pending planning approval

Rows / page: 15 30 50

< Previous 1 2 Next >

A normal click on a patient row will display the full details of the case

Which patients are displayed in the list ?

Symbios Connect Dr CHARNLEY John ▾

Surgical cases

 New case 

Date quick filter

July 2022						
Mo	Tu	We	Th	Fr	Sa	Su
27	28	29	30	1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31
1	2	3	4	5	6	7

Today

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN®	Shipped
	HARI Mata	KNEE	R	ORIGIN®	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN®	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN®	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available



- **Surgeon account:** visibility only on the **surgeon's patients**
- **Sales rep account:** visibility on **all patients that are operated in hospitals** that are linked to the sales rep's account
- **Secretary account:** visibility on **all patients that are operated by surgeons** that are linked to the secretary's account

Quick filter by status – pending approval cases



Deactivate all other status filters and activate “Pending approvals” filter

Pr CHARNLEY John ▾

Symbios Connect

Surgical cases

New case

Search



Date quick filter

« February 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	1	2	3	4	5	6
7	8	9	10	11	12	13

Today

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
06.09.2022	MONROE Marilyn	KNEE	R	ORIGIN®	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN®	Pending planning approval
15.09.2022	CHURCHILL Winston	HIP	R	INDIVIDUAL HIP™	Pending planning approval
04.10.2022	KELLY Grace	KNEE	L	ORIGIN®	Pending planning approval

< Previous 1 Next >

Status quick filters

Active cases [19]

Pending approvals [4]

In process [14]

Delayed [1]

Inactive cases [1]

Operated [1]

Cancelled [0]

Quick filter by status – Delayed cases

**Deactivate all other status filters
And activate “Delayed” filter**

Pr CHARNLEY John ▾

Symbios Connect

Surgical cases

New case

Search

Date quick filter

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HARI Mata	KNEE	R	ORIGIN®	■ Waiting for CT scan

HARI Mata KNEE RIGHT (126361) ✕

Related order

ORIGIN®

Custom knee prosthesis

Case ID 126361

Planned implants

- ORIGIN CR Femur Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented

Ordered products

- ORIGIN CR Femur Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented
- ORIGIN KNEE-PLAN Femoral Cut Guide
- ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

Quick filter by status – Operated cases



Make sure to select a date in the date filter that is prior to the surgery date

Sym

Pr CHARNLEY John ▾

Surgical cases

New case



Deactivate all other status filters and activate “Operated” filter



Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.02.2022	TAYLOR Elizabeth	KNEE	R	ORIGIN®	Operated

TAYLOR Elizabeth
KNEE RIGHT
(126377) ×

Related order

ORIGIN® Modular

Custom knee prosthesis

Case ID 126377

Planned implants

- ORIGIN PS Femur Cemented
- ORIGIN PS Fixed Tibia Modular Cemented
- ORIGIN PS Fixed Insert
- ORIGIN Patella Cemented

Ordered products

- ORIGIN PS Femur Cemented
- ORIGIN PS Fixed Tibia Modular Cemented
- ORIGIN PS Fixed Insert
- ORIGIN Patella Cemented
- ORIGIN KNEE-PLAN Femoral Cut Guide
- ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

Status quick filters

Active cases (19)

- Pending approvals (4)
- In process (14)
- Delayed (1)

Inactive cases (1)

- Operated (1)
- Cancelled (0)
- Suspended (0)

Filter by type of orders/product



Click on the “Related order” filter and select which product type to filter

Symbios Connect Pr CHARNEY John ▾

Surgical cases

New case Search

Date quick filter

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	FORD Henri	KNEE	R	All	Shipped
	HARI Mata	KNEE	R	ORIGIN®	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	HIP-PLAN	Manufacturing in progress
	BAKER Joséphine	KNEE	L	INDIVIDUAL HIP™	Ready for delivery
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN®	Planning approved
	MONROE Marilyn	KNEE	R	ORIGIN®	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN®	Pending planning approval
04.10.2022	KELLY Grace	KNEE	L	ORIGIN®	Pending planning approval
11.10.2022	DOE John	KNEE	L	ORIGIN®	Waiting for CT scan
	WASHINGTON George	KNEE	R	ORIGIN®	Waiting for CT scan



Click on ‘All’ to reset a column filter and display all results

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]

Filter by status (exhaustive list of status)



Click on the "Status" filter and select which status to filter

Pr CHARNEY John ▾

Symbios Connect

Surgical cases

New case

Search



Date quick filter

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	All
	FORD Henri	KNEE	R	ORIGIN®	Shipped

< Previous 1 Next >

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]

- All
- Shipped
- CT scan available
- CT scan validated
- Manufacturing in progress
- Pending planning approval
- Planning approved

Make sure to reset filters to display normally

Symbios Connect

Pr CHARNLEY John ▾

Surgical cases

New case

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	All
	FORD Henri				Shipped
					CT scan available
					CT scan validated

Status ▾

Date quick filter

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]

Warning: If display is not 'normal' make sure to check for potential activated filter(s). Check for :

1. Is any 'quick status filter' activated ?
2. Is any column filter activated ? Check for the icon above on top of a column to verify this.
3. Is the date filter preventing some case to be displayed ?

Search a case



Type in the search textbox to search a case. The results appear in real-time in a window below the textbox.

Symbios Connect

Surgical cases

Search criteria: Name of patient, case ID

New case

Date quick filter

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]

Surgery date	Patient	Joint				
23.08.2022	HEPBURN Audrey	HIP				
	FORD Henri	KNEE				
	HARI Mata	KNEE	R	ORIGIN®	<input type="checkbox"/>	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN®	<input checked="" type="checkbox"/>	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN®	<input checked="" type="checkbox"/>	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	<input checked="" type="checkbox"/>	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	<input checked="" type="checkbox"/>	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN®	<input checked="" type="checkbox"/>	Planning approved
	CHANEL Coco	HIP	L	INDIVIDUAL HIP™	<input checked="" type="checkbox"/>	Planning in progress
	MONROE Marilyn	KNEE	R	ORIGIN®	<input type="checkbox"/>	Pending planning approval

< Previous 1 2 Next >

Click on a case in the list to display its details

DOE John

LEFT KNEE

[126378]

Waiting for CT scan

Patient information

Administrative information

Patient name	DOE John
Date of birth	07.02.1968
Gender	Male
Comments	

Medical information

Joint to operate	Knee
Side	Left
Case ID	126378
Preop. extension	5° Flessum
Preop. flexion	30°
Preop. deformity	Varus 12°
Passive correction	Total

[← Back to surgical cases](#)

[✎ Edit](#)

[⌛ Suspend surgery](#)

[⌛ Cancel surgery](#)

Products to order

ORIGIN®

Planned implants	ORIGIN CR Femur Cemented ORIGIN CR Fixed Tibia Monobloc Cemented ORIGIN CR Fixed Insert ORIGIN Patella Cemented
Ordered products	ORIGIN CR Femur Cemented ORIGIN CR Fixed Tibia Monobloc Cemented ORIGIN CR Fixed Insert ORIGIN Patella Cemented ORIGIN KNEE-PLAN Femoral Cut Guide ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

[✎ Edit](#)

CT scan exam and surgery date

Radiology center	Wigan Albert Edward Infirmary - WN1 2NN Wigan
CT exam date	09.08.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

[✎ Edit](#)

Create a new case

How to place a custom product order

Click on “new case” button

Symbios Connect Dr CHARNLEY John ▾

Click on « new case » button

Surgical cases

New case

Q

Date quick filter	Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾																																																	
<div style="border: 1px solid #ccc; padding: 5px;"> <p>« July 2022 »</p> <table style="width: 100%; text-align: center; font-size: 10px;"> <tr> <td>Mo</td><td>Tu</td><td>We</td><td>Th</td><td>Fr</td><td>Sa</td><td>Su</td> </tr> <tr> <td>27</td><td>28</td><td>29</td><td>30</td><td>1</td><td>2</td><td>3</td> </tr> <tr> <td>4</td><td>5</td><td>6</td><td>7</td><td>8</td><td>9</td><td>10</td> </tr> <tr> <td>11</td><td>12</td><td>13</td><td>14</td><td>15</td><td>16</td><td>17</td> </tr> <tr> <td>18</td><td>19</td><td>20</td><td>21</td><td>22</td><td>23</td><td>24</td> </tr> <tr> <td>25</td><td>26</td><td>27</td><td style="background-color: #e0f2f1;">28</td><td>29</td><td>30</td><td>31</td> </tr> <tr> <td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td> </tr> </table> <p style="text-align: right; font-size: 10px;">Today</p> </div>	Mo	Tu	We	Th	Fr	Sa	Su	27	28	29	30	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	1	2	3	4	5	6	7	23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	Mo	Tu	We	Th	Fr	Sa	Su																																																
	27	28	29	30	1	2	3																																																
	4	5	6	7	8	9	10																																																
	11	12	13	14	15	16	17																																																
	18	19	20	21	22	23	24																																																
	25	26	27	28	29	30	31																																																
	1	2	3	4	5	6	7																																																
			FORD Henri	KNEE	R	ORIGIN®	Shipped																																																
			HARI Mata	KNEE	R	ORIGIN®	Waiting for CT scan																																																
	30.08.2022	DESCARTES René	KNEE	R	ORIGIN®	Manufacturing in progress																																																	
		BAKER Joséphine	KNEE	L	ORIGIN®	Ready for delivery																																																	
		CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available																																																	
		ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available																																																	
	06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN®	Planning approved																																																	
		CHANEL Coco	HIP	L	INDIVIDUAL HIP™	Planning in progress																																																	
		MONROE Marilyn	KNEE	R	ORIGIN®	Pending planning approval																																																	
	13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN®	Pending planning approval																																																	
		MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available																																																	
	15.09.2022	CHURCHILL Winston	HIP	R	INDIVIDUAL HIP™	Pending planning approval																																																	
		FITZGERALD Ela	HIP	R	HIP-PLAN	Waiting for CT scan																																																	
	04.10.2022	KELLY Grace	KNEE	L	ORIGIN®	Pending planning approval																																																	

Rows / page: 15 30 50

< Previous 1 2 Next >

Status quick filters

Active cases (18)

- Pending approvals (4)
- In process (13)
- Delayed (1)

Inactive cases (1)

- Operated (1)
- Cancelled (0)
- Suspended (0)
- To complete (0)

Enter patient information

Dr CHARNLEY John ▾

Symbios Connect

New surgical case

1 Patient information

Administrative information

Gender * Male Female

Medical information

Joint to operate * Hip Knee

* Mandatory

[Back to surgical cases](#)[View orders lead time...](#)[Save as incomplete](#)[Continue](#)[2 Products to order](#)[3 CT scan exam and surgery date](#)[4 Confirm case](#)

Enter patient information

Dr CHARLEY John ▾

Symbios Connect


Fill in patient's information

New surgical case

1 Patient information

Administrative information

DOE

John

07 Feb 1968

Gender * Male Female

Medical information

Joint to operate * Hip Knee

RIGHT KNEE  LEFT KNEE

* Mandatory

Comments

Preop. extension
 15° 30°
 Recurvatum Flessum

Preop. flexion
 30° 140°

Preop. deformity (loaded conditions)
 30° 30°
 Valgus Varus

Passive correction
 Total Partial Fixed

Preoperative clinical information

New surgical case

1 Patient information

Administrative information


Gender * Male Female

Preoperative clinical info

Highly recommended that surgeon provides them, as it significantly improves the quality of planning

* Mandatory


Preop. extension



Recurvatum

Flessum

Preop. flexion



Preop. deformity (loaded conditions)



Valgus

Varus

Passive correction



[Back to surgical cases](#)

[View orders lead time...](#)

[Save as incomplete](#)

[Continue](#)

2 Products to order

3 CT scan exam and surgery date

Real-time product lead time

New surgical case

1 Patient information



Click to view delivery time for custom orders

Administrative information

DOE

John

07 Feb 1968

Gender * Male Female

Medical information

Joint to operate * Hip Knee

RIGHT KNEE



LEFT KNEE

* Mandatory

Comments

Preop. extension

15° 5° 30°
Recurvatum Flessum

Preop. flexion

30° 140°

Preop. deformity (loaded conditions)

30° 12° 30°
Valgus Varus

Passive correction

N/D Total Partial Fixed

Back to surgical cases

View orders lead time...

Save as incomplete

Continue

2 Products to order

3 CT scan exam and surgery date

Real-time product lead time

1

Select the custom-made product for which you would like to check the delivery time

1 Patient info

Administr

DOE

John

07 Feb 196

Gender *


Medical i

Joint to o



INDIVIDUAL HIP™

Personalized hip femoral stem



KNEE-PLAN®

Preoperative 3D planning and custom instruments



ORIGIN®

Custom knee prosthesis

2

Select date of CT scan transfer to Symbios

an exam

today << < August 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

* Mandator

Back to surgical cases

2 Products to order

3 CT scan exam and su

3

View possible dates for surgery

Possible surgery date
on CT exam date

today << < October 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

s incomplete

Continue

Real-time product lead time



Light grey dates could be available for surgery if CT was sent at earlier date than the one currently selected

If you select a light grey date, the left calendar will propose the latest dates for CT to be transferred to Symbios

today

Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

OK

Save as incomplete (optional)



Click here to save patient case and complete it later

New surgical case

1 Patient information

Administrative information

DOE

John

07 Feb 1968

Gender * Male Female

Medical information

Joint to operate * Hip Knee

RIGHT KNEE



LEFT KNEE

* Mandatory

[Back to surgical cases](#)

[View orders lead time...](#)

2 Products to order

3 CT scan exam and surgery date

Comments

Preop. extension

15° 5° 30°

Recurvatum Flessum

Preop. flexion

30° 140°

Preop. deformity (loaded conditions)

30° 12° 30°

Valgus Varus

Passive correction

N/D Total Partial Fixed

Save as incomplete

Continue

How to find incomplete orders in patient list ?


Symbios Connect Dr CHARNLEY John ▾

 **Incomplete order**

Surgical cases Search 

Activate the filter
« to complete »

Deactivate all other filters to easily find the incomplete order(s)

Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
MONROE Marilyn	KNEE	R	ORIGIN@	 Pending planning approval
EINSTEIN Albert	KNEE	R	ORIGIN@	 Pending planning approval
MANDELA Nelson	HIP	L	HIP-PLAN	 CT scan available
CHURCHILL Winston	HIP	R	INDIVIDUAL HIP™	 Pending planning approval
FITZGERALD Ela	HIP	R	HIP-PLAN	 Waiting for CT scan
KELLY Grace	KNEE	L	ORIGIN@	 Pending planning approval
WAYNE John	HIP	L	INDIVIDUAL HIP™	 CT scan validated
11.10.2022 VAN BEETHOVEN Ludwig	HIP	L	INDIVIDUAL HIP™	 Waiting for CT scan
WASHINGTON George	KNEE	R	ORIGIN@	 Waiting for CT scan
- DOE John	KNEE	L		 Pending order

Inactive cases (2)

Operated (1)

Cancelled (0)

Suspended (0)

To complete (1)

< Previous 1 Next >

Rows / page: 15 30 **50**

How to create an ORIGIN® order ?



Click on the **ORIGIN®** button

Administrative information

Patient name DOE John
Date of birth 07.02.1968

Medical information

Joint to operate	Knee
Side	Left
Preop. extension	5° Flessum
Preop. flexion	30°
Preop. deformity	Varus 12°
Passive correction	Total

[Edit](#)

2 Products to order



KNEE-PLAN®

Preoperative 3D planning and custom instruments



ORIGIN®

Custom knee prosthesis

[Back to surgical cases](#)

[View orders lead time...](#)

[Continue](#)

3 CT scan exam and surgery date

4 Confirm case

Select ORIGIN[®] components and options

1

Select product options and components

Preoperative CT planning and custom instruments

Please select the implants you want to plan for this patient:



ORIGIN CR Femur Cemented
Custom-made cemented femoral component

✓ Change

ORIGIN CR Femur Cemented

ORIGIN PS Femur Cemented



ORIGIN CR Fixed Tibia Monobloc Cemented
Custom-made cemented monobloc tibial component for fixed platform

Change ▾



ORIGIN CR Fixed Insert
Custom-made tibial fixed insert



Please select if you want to plan a patella for this patient:



ORIGIN Patella Cemented
Cemented patellar component with dome. Polyethylene (UHMWPE-ISO 5

2

Select custom-guide design according to the chosen surgical approach

Please select the custom cutting-guides that you want to order for this patient:



ORIGIN KNEE-PLAN Femoral Cut Guide
Custom-made femoral cut guide. Polyamide PA2200.
Delivered sterile with femoral bone model



ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD
Custom-made tibial cut guide for internal approaches. Polyamide PA2200.
Delivered sterile with tibial bone model

Surgical approach

Medial (default)

Lateral

Recommended to adapt surgeon's preferences for products he regularly uses in his settings (see section 8). The preferences will be automatically loaded at this step

How to create an INDIVIDUAL HIP[®] stem order ?



Click on the INDIVIDUAL HIP[®] button

Comments

Surgical information

Joint to operate


Side

Hip

Right

Edit

2 Products to order



[Back to surgical cases](#)

[View orders lead time...](#)

[Continue](#)

3 CT scan exam and surgery date

4 Confirm case

Select INDIVIDUAL HIP® stem/rasp design options

1

Select stem options



INDIVIDUAL HIP TiHA

Personalized cementless stem. Titanium alloy (Ti6Al4V-ISO 5832-3), coated with porous titanium and hydroxyapatite. 12/14 5°40' taper.

Change ▾

Design options

Coating

- Partial HA coating
- Full HA coating

Collar

- No collar
- Normal collar



INDIVIDUAL HIP Rasp

Personalized piqued rasp. Stainless steel (316L - NF S 94-090). 12/14 5°40' taper.

Change ▾

Design options

Impaction hole

- With thread (use with INDIVIDUAL HIP™ rasp handle)
- Without thread (use with INDIVIDUAL HIP™ stem impactor)



BIOLOX Delta Head

Ceramic head (Al2O3 + ZrO2-ISO 6474-2), compatible with 12/14 5°40' taper.

Change ▾



APRIL Ceramic

Cementless pressfit acetabular cup. Titanium alloy (Ti6Al4V-ISO 5832-3). Coatings: Porous titanium and hydroxyapatite.

Change ▾



BIOLOX Delta Insert

Truncated conical insert. Ceramic (Al2O3 + ZrO2-ISO 6474-2).

2

Select rasp options

Recommended to adapt surgeon's preferences for products he regularly uses in his settings (see section 8). The preferences will be automatically loaded at this step

How to create a HIP-PLAN® case ?



Click on the HIP-PLAN® button

Comments

Medical information

Joint to operate



Side

Hip

Right

 Edit

2 Products to order



[Back to surgical cases](#)

[View orders lead time...](#)

[Continue](#)

3 CT scan exam and surgery date

4 Confirm case

Select standard hip implants to plan

1

Select femoral components



HIP-PLAN®

3D preoperative planning



INDIVIDUAL HIP™

Personalized hip femoral stem



Recommended to adapt surgeon's preferences for products he regularly uses in his settings (see section 8). The preferences will be automatically loaded at this step

2

Select acetabular components

Please select the femoral components you want to plan for this patient:



SPS Evolution

Cementless anatomic femoral stem. Titanium alloy (Ti6Al4V-ISO 5832-3). Coatings: Porous titanium and hydroxyapatite. 12/14 5°40' taper. CCD shaft angle 129°.

Change ▾



BIOLOX Delta Head

Ceramic head (Al₂O₃ + ZrO₂-ISO 6474-2), compatible with 12/14 5°40' taper.

Change ▾

Please select the acetabular components you want to plan for this patient:



APRIL Ceramic

Cementless pressfit acetabular cup. Titanium alloy (Ti6Al4V-ISO 5832-3). Coatings: Porous titanium and hydroxyapatite.

Change ▾



BIOLOX Delta Insert


Truncated conical insert. Ceramic (Al₂O₃ + ZrO₂-ISO 6474-2).

How to create a KNEE-PLAN® case ?



Click on the KNEE-PLAN® button

2 Products to order



Joint to operate	Knee
Side	Left
Preop. extension	5° Flexum
Preop. flexion	30°
Preop. deformity	Varus 12°
Passive correction	Total

[Back to surgical cases](#)

[View orders lead time...](#)

[Continue](#)

3 CT scan exam and surgery date

4 Confirm case

Select FIRST implants and KNEE-PLAN® options

1

Select product options and components

Preoperative

Please select the implants you want to plan for this patient:



FIRST Femur Cemented

Cemented femoral component. Cobalt-Chrome (CoCrMo-ISO 5832-4).

Change ▾



FIRST Fixed Tibia Monobloc Cemented

Cemented monobloc tibial component for fixed platform. Titanium alloy (Ti6Al4V-ISO 5832-3).

Change ▾



FIRST Fixed Insert

Tibial fixed insert. Polyethylene (UHMWPE-ISO 5834-2).

Please select if you want to plan a patella for this patient:



FIRST Patella

Cemented patella component. Polyethylene (UHMWPE-ISO 5834-2).

Change ▾

2

Select custom-guide design according to the chosen surgical approach

Please select the custom cutting-guides that you want to order for this patient:



FIRST KNEE-PLAN Femoral Cut Guide - STD

Custom-made femoral cut guide. Compatible with STD 4-in-1 femoral cuts guide. Polyamide PA2200.

Delivered sterile with femoral bone model



FIRST KNEE-PLAN Tibial Cut Guide - INT STD

Custom-made tibial cut guide for internal approaches. Polyamide PA2200.

Delivered sterile with tibial bone model

Surgical approach

 Medial (default)

 Lateral

FIRST KNEE-PLAN tibial cut guides' design is optimized to fit the surgical approach you plan to use



Recommended to adapt surgeon's preferences for products he regularly uses in his settings (see section 8). The preferences will be automatically loaded at this step

Select medical imaging centre for the CT exam

3 CT scan exam and surgery date

Select medical imaging centre for CT scan exam



Click here to search and select a medical imaging centre

WRIGHTINGTON HOSPITAL

Appley Bridge

Theatre Stores
WN6 9EP Appley Bridge
UNITED KINGDOM

Tel: +44 1942 244000

Fax:

Select date for CT scan exam

today << < August 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

Select date of surgical procedure

today << < October 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

[Back to surgical cases](#)

[Continue](#)

Select medical imaging centre for the CT exam

1 Fill in search parameter(s) and click on magnifier

Search for a radiology centre

Enter the country and locality information to refine the search for a radiology centre that is a Symbios Connect partner

Name of centre

United Kingdom

London



Radiology center

Newham Hospital

Glen Road
E13 8SL London
UNITED KINGDOM

Email:
Tel: 020 7476 4000
Fax:

2 Select the medical imaging centre

The Whittington CT

Magdala Avenue
N19 5NF London
UNITED KINGDOM

Email:
Tel:
Fax:

London Bridge Hospital

Tooley Street 27
SE12PR LONDON
UNITED KINGDOM

Email:
Tel: 20 72 34 27 73
Fax:

Guys & St Thomas' Nhs Trust

Email: stt.nhs.uk

3

Click here to select



Click here to add as favourite

Cancel

Add to favourites

Select

Back to surgical cases

Continue

Select hospital for surgery

3 CT scan exam and surgery date

Select medical imaging centre for CT scan exam

Select hospital for surgery

WRIGHTINGTON HOSPITAL

Appley Bridge

Theatre Stores
WN6 9EP Appley Bridge
UNITED KINGDOM

Tel: +44 1942 244000
Fax:

Click here to change the hospital (in case surgeon operates in several hospitals)

Select date of surgical procedure

today << < August 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

today << < October 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

Select dates for CT and surgery

3 CT scan exam and surgery date

1

Select the date on which the CT scan will be transferred to Symbios

WN1 2NN Wigan
UNITED KINGDOM

Fax:

Favourites

Select date for CT scan exam

today << < August 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

2

Select the date on which the surgery should be planned

Information on first available dates for the surgery will be provided in this calendar in white

Select ho

WRIGHTIN
Appley Br
Theatre St
WN6 9EP Ap
UNITED KIN

Select date of surgical procedure

today << < October 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

How to use real-time delivery information



It is not possible to select dates that are indicated in black, as Symbios cannot deliver for the surgery if the CT is provided on the date selected on the left



Light grey dates could be available for surgery if CT was sent at earlier date than the one currently selected

Select date for CT scan exam

today << < August 2022 > >>						
Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

Select date of surgical procedure

today << < October 2022 > >>						
Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

How to use real-time delivery information



If you select a light grey date, the left calendar will propose **the latest date for CT to be transferred to Symbios, which is the last date in white** before the grey dates

Select date for CT scan exam

today	<<	<	August 2022	>	>>	
Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

Select date of surgical procedure

today	<<	<	October 2022	>	>>	
Sun	Mon	Tue	Wed	Thu	Fri	Sat
			5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

[Back to surgical cases](#)

[Continue](#)

Finalize the order

Gender Male

Comments

Preop. extension 5° Flessum
 Preop. flexion 30°
 Preop. deformity Varus 12°
 Passive correction Total

[Edit](#)

2 Products to order



ORIGIN®

Planned implants
 ORIGIN CR Femur Cemented
 ORIGIN CR Fixed Tibia Monobloc Cemented
 ORIGIN CR Fixed Insert
 ORIGIN Patella Cemented

Ordered products
 ORIGIN CR Femur Cemented
 ORIGIN CR Fixed Tibia Monobloc Cemented
 ORIGIN CR Fixed Insert
 ORIGIN Patella Cemented
 ORIGIN KNEE-PLAN Femoral Cut Guide
 ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD



Click on “confirm order” to finalize

3 CT scan exam and surgery date

Radiology center Wigan Albert Edward Infirmary - WN1 2NN Wigan
 CT exam date 09.08.2022

Hospital WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
 Surgery date 11.10.2022

[Edit](#)

4 Confirm case

By clicking on **Confirm case**, I confirm that the above information is correct, and confirm that I have read and accepted Symbios Connect terms of use and conditions.

[Back to surgical cases](#)

[Confirm order](#)

New surgical case registered

Gender Male

New surgical case

Your order has been created with ID:126378

[Back to the order list](#)

[Order details](#)

[Edit](#)

2 Products to order



ORIGIN®

Ordered products

ORIGIN CR Fixed insert
ORIGIN Patella Cemented
ORIGIN CR Femur Cemented
ORIGIN CR Fixed Tibia Monobloc Cemented
ORIGIN CR Fixed Insert
ORIGIN Patella Cemented
ORIGIN KNEE-PLAN Femoral Cut Guide
ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

3 CT scan exam and surgery date

Radiology center Wigan Albert Edward Infirmary - WN1 2NN Wigan
CT exam date 09.08.2022

[Edit](#)

Hospital WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date 11.10.2022

4 Confirm case

By clicking on **Confirm case**, I confirm that the above information is correct, and confirm that I have read and accepted Symbios Connect terms of use and conditions.

[Back to surgical cases](#)

[Confirm order](#)

Planning approval

How to approve the preop planning for custom products

Method n°1 – Approve via the email



An email is sent to the surgeon when the planning is completed

Planning pending ORIGIN® approval - EINSTEIN Albert - (126369)



• Symbios Connect <connect@symbios.ch>

À : sc.validation1@symbios.ch

Hier à 16:49



Télécharger tout · Prévisualiser tout

1

Control planning with the attached PDF

Hello Pr CHARNEY,



The planning report for this ORIGIN® case is now available and you can view it by opening the email attachment, or by logging on to [Symbios Connect](#).

Our design teams are awaiting your approval of the proposed planning, before launching production of the custom-made products associated with this case.



Approval deadline:
16.08.2022

[Approve the planning directly](#)

[View case in Symbios Connect](#)

Reminder: completed and updated information for this case is provided in detail below.

Patient information

Administrative information

Patient name EINSTEIN Albert
Date of birth 18.04.1955
Gender Male
Case ID **126369**
Surgeon Pr CHARNEY John
Comments

Clinical information

Joint **Knee**
Side **Right**

Products to order



Planned implants ORIGIN CR Femur Cemented
ORIGIN CR Fixed Tibia Monobloc Cemented
ORIGIN CR Fixed Insert
ORIGIN Patella Cemented

Ordered products ORIGIN CR Femur Cemented
ORIGIN CR Fixed Tibia Monobloc Cemented

Method n°1 – Approve via the email

2

Click on the **green link** to approve the planning

OR

2

Click on the **grey link** if want to ask a question or modify

Planning pending ORIGIN® approval - EINSTEIN Albert - (126369)



Symbios Connect <connect@symbios.ch>

À: sc.validation1@symbios.ch



Télécharger tout · Prévisualiser tout

Hello Pr CHARNEY,

The planning report for this ORIGIN® case is now available and you can view it by opening the email attachment, or by logging on to [Symbios Connect](#).

Our design teams are awaiting your approval of the proposed planning, before launching production of the custom-made products associated with this case.



Approval deadline:
16.08.2022

[Approve the planning directly](#)

[View case in Symbios Connect](#)

Reminder: completed and updated information for this case is provided in detail below.

Patient information

Administrative information

Patient name EINSTEIN Albert
Date of birth 18.04.1955

Clinical information

Joint **Knee**
Side **Right**

Make sure to approve before the approval deadline

Products to order



Planned implants
ORIGIN CR Femur Cemented
ORIGIN CR Fixed Tibia Monobloc Cemented
ORIGIN CR Fixed Insert
ORIGIN Patella Cemented

Ordered products
ORIGIN CR Femur Cemented
ORIGIN CR Fixed Tibia Monobloc Cemented



Method n°2 – Approve via Symbios Connect


EINSTEIN Albert

RIGHT KNEE
[126369]

Pending planning approval

Patient information

Administrative information	Medical information
Patient name	EINSTEIN Albert
Date of birth	18.04.1955
Gender	Male
	Joint to operate
	Side
	Case ID
	Knee
	Right
	126369

[← Back to surgical cases](#)

[⏸ Suspend surgery](#)

[⏹ Cancel surgery](#)

Products to order



ORIGIN®

**Click on “Approve planning”
to approve formally the case**

Ordered products

- ORIGIN Patella Cemented
- ORIGIN CR Femur Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented
- ORIGIN KNEE-PLAN Femoral Cut Guide
- ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

⌚ 2 days
remaining to
approve planning

[📄 Download latest
planning.](#)

[✅ Approve planning.](#)

[⚙ Modify planning.](#)

[❓ Ask a question to
engineer.](#)

[🕒 View planning history.](#)

Make sure to approve before the approval deadline

[✎ Edit](#)

Hospital

WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

Ask a question to the engineer

 EINSTEIN Albert RIGHT KNEE [126369]

Pending planning approval

Patient information

[← Back to surgical cases](#)

Administrative information

Patient name EINSTEIN Albert
Date of birth 18.04.1955
Gender Male

Medical information

Joint to operate Knee
Side Right
Case ID 126369

 [Suspend surgery](#)
 [Cancel surgery](#)

Products to order



**Click on “ask a question”
to raise question to the engineer
who planned the case**



ORIGIN®

products ORIGIN CR Fixed Tibia Monobloc Cemented
ORIGIN CR Fixed Insert
ORIGIN Patella Cemented
ORIGIN KNEE-PLAN Femoral Cut Guide
ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

 **2 days**
remaining to
approve planning
 [Download latest
planning.](#)
 [Approve planning.](#)
 [Modify planning.](#)
 [Ask a question to
engineer.](#)
 [View planning history.](#)

CT scan exam and surgery date

Radiology center Radiologie Stuttgart - 70173 Stuttgart
CT exam date 01.07.2022

Hospital WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

 [Edit](#)

Ask a question to the engineer

The screenshot shows a web interface for planning approval. A modal window titled "Ask a question to engineer" is open, containing a text input field and "Cancel" and "Submit" buttons. A green callout box with a lightbulb icon points to the input field, containing the text "Fill in question to engineer and click 'submit'".

Ask a question to engineer

Please type here any question you want to ask the Custom engineer. We will take this request into consideration and reply as soon as possible.

Enter your question here

Cancel Submit

Fill in question to engineer and click "submit"

2 days remaining to approve planning

- Download latest planning
- Approve planning
- Modify planning
- Ask a question to engineer
- View planning history

[Edit](#)

Patient information

Administrative

Products to order

ORIGIN®

Ordered products

- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented
- ORIGIN CR Femur Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Insert

CT scan exam and surgery date

Radiology center	Radiologie Stuttgart - 70173 Stuttgart
CT exam date	01.07.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appleby Bridge
Surgery date	13.09.2022

Modify planning

 EINSTEIN Albert RIGHT KNEE [126369]

Pending planning approval

Patient information

[← Back to surgical cases](#)

Administrative information

Patient name EINSTEIN Albert
Date of birth 18.04.1955
Gender Male

Medical information

Joint to operate Knee
Side Right
Case ID 126369

 [Suspend surgery](#)
 [Cancel surgery](#)

Comments



**Click on “Modify planning”
to request a modification of the
preoperative planning / design**

Products to order



ORIGIN®

Planned implants
ORIGIN CR
ORIGIN CR
ORIGIN CR
ORIGIN Pat

Ordered products
ORIGIN CR
ORIGIN CR Fixed Tibia Monobloc Cemented
ORIGIN CR Fixed Insert
ORIGIN Patella Cemented
ORIGIN KNEE-PLAN Femoral Cut Guide
ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

 **2 days**
remaining to
approve planning
 [Download latest
planning.](#)
 [Approve planning.](#)
 [Modify planning.](#)
 [Ask a question to
engineer.](#)
 [View planning history.](#)

CT scan exam and surgery date

Radiology center Radiologie Stuttgart - 70173 Stuttgart
CT exam date 01.07.2022

Hospital WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

 [Edit](#)

Modify planning

The screenshot shows a web application interface for planning approval. A modal dialog box titled "Modification request" is open, prompting the user to enter modifications. The dialog contains a text input field labeled "Modifications requested" and "Submit" and "Cancel" buttons. A green arrow points from a callout box to the input field. The background page shows patient information, a list of ordered products (including ORIGIN CR Fixed Insert, Patella Cemented, Femur Cemented, and Tibia Monobloc Cemented), and surgery dates. A callout box with a lightbulb icon contains the instruction: "Fill in modification request to the engineer and click 'submit'".

Modification request

Please type here the modifications to the planning that you would request from our custom engineers. We will take this request into consideration and make another proposal to you as soon as possible.

Modifications requested

Cancel Submit

Fill in modification request to the engineer and click "submit"

2 days remaining to approve planning

- [Download latest planning.](#)
- [Approve planning.](#)
- [Modify planning.](#)
- [Ask a question to engineer.](#)
- [View planning history.](#)

[Back to surgical cases](#)

[Suspend surgery.](#)

[Cancel surgery.](#)

[Edit](#)

Patient information

Administrative

Products to order

ORIGIN®

Ordered products

- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented
- ORIGIN CR Femur Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented

CT scan exam and surgery date

Radiology center	Radiologie Stuttgart - 70173 Stuttgart
CT exam date	01.07.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date	13.09.2022

Modify a case

Display the case details to modify some elements

DOE John

LEFT KNEE

[126378]

Waiting for CT scan

Patient information

Administrative information

Patient name	DOE John
Date of birth	07.02.1968
Gender	Male

Comments

**Click on "Edit" button(s)
to modify case information**

Medical information

Joint to operate	Knee
Side	Left
Case ID	126378
Preop. extension	5° Flexum
Preop. flexion	30°
Preop. deformity	Varus 12°
Passive correction	Total

[← Back to surgical cases](#)

[✎ Edit](#)

[⏸ Suspend surgery](#)

[⏹ Cancel surgery](#)

Products to

ORIGIN®

Planned implants	ORIGIN CR Femur Cemented ORIGIN CR Fixed Tibia Monobloc Cemented ORIGIN CR Fixed Insert ORIGIN Patella Cemented
Ordered products	ORIGIN CR Femur Cemented ORIGIN CR Fixed Tibia Monobloc Cemented ORIGIN CR Fixed Insert ORIGIN Patella Cemented ORIGIN KNEE-PLAN Femoral Cut Guide ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

✎ [Edit](#)

CT scan exam and surgery date

Radiology center	Wigan Albert Edward Infirmary - WN1 2NN Wigan
CT exam date	09.08.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

✎ [Edit](#)

Edit patient information

Patient information

Administrative information

 Gender * Male Female



Modify and confirm by clicking on "Confirm" button

 RIGHT KNEE LEFT KNEE


Extension

Preop. flexion

Flexum

Preop. deformity (loaded conditions)

Valgus Varus

Passive correction

Partial Fixed

* Mandatory

Products to order



Planned implants

- ORIGIN CR Femur Cemented
- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented

Edit product information



Custom-made cemented femoral component Change ▾



ORIGIN CR Fixed Tibia Monobloc Cemented
Custom-made cemented monobloc tibial component for fixed platform Change ▾



Custom-made femoral cut guide. Polyamide PA2200.
Delivered sterile with femoral bone model



ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD
Custom-made tibial cut guide for internal approaches. Polyamide PA2200.
Delivered sterile with tibial bone model

Medial (default)
 Lateral

Modify and confirm by clicking on “Confirm” button

Please select if you want to plan a patella for this patient:



ORIGIN Patella Cemented
Cemented patellar component with medialized dome. Polyethylene [UHMWPE-ISO 5834-2].

CT scan exam and surgery date

Cancel Confirm

Radiology center Wigan Albert Edward Infirmary - WN1 2NN Wigan
 CT exam date 09.08.2022

Edit CT/surgery information

CT scan exam and surgery date

Select medical imaging centre for CT scan exam

Wigan Albert Edward Infirmary

Wigan

Wigan Lane

Wigan

UNITED KINGDOM

Tel: +44 1942 244000

Fax:

Select hospital for surgery

WRIGHTINGTON HOSPITAL

Appley Bridge

Theatre Stores

WN6 9EP Appley Bridge

UNITED KINGDOM

Tel: +44 1942 244000

Fax:

Change hospital



Modify and confirm by clicking on "Confirm" button

Select

today << < August 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
31	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

Select date of surgical procedure

today << < October 2022 > >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
				6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

Cancel **Confirm**

Not always possible to edit – dependence on status



Not possible to edit patient or product information as soon as the design process has started

Date of birth	24.01.1965	Side	Right
Gender	Male	Case ID	126371
Comments		Weight	96 kg
		Size	171 cm

Pending planning approval

← [Back to surgical cases](#)

-  [Suspend surgery](#)
-  [Cancel surgery](#)

Products to order

Planned implants	INDIVIDUAL HIP Cemented Normal collar HILOCK LINE INLOCK X
Ordered products	INDIVIDUAL HIP Cemented Normal collar INDIVIDUAL HIP Compactor NT Threaded impaction hole





Still possible to modify the surgery date (with limitations)


 **6 days**
remaining to approve planning

-  [Download latest planning](#)
-  [Approve planning](#)
-  [Modify planning](#)
-  [Ask a question to engineer](#)
-  [View planning history](#)

CT scan exam and surgery date

Radiology center	King's Mill Hospital - NG17 4JL Sutton-In-Ashfield
CT exam date	01.07.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

✎ Edit



Modify the date of surgery



can exam and surgery date

Possible to modify the surgery date

Limitations usually are linked to the **minimum lead-time** necessary for the production of product. Always possible to postpone the date.

Sun	Mon	Tue	Wed	Thu	Fri	Sat
26	27	28	29	30	1	2
3	4	5	6	7	8	9
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31	1	2	3	4	5	6

Select hospital for surgery

WRIGHTINGTON HOSPITAL

Appley Bridge

Theatre Stores
WNG 9EP Appley Bridge
UNITED KINGDOM

Tel: +44 1942 244000

Fax:

Change hospital

Select date of surgical procedure

today << September 2022 >>

Sun	Mon	Tue	Wed	Thu	Fri	Sat
	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	1
2	3	4	5	6	7	8

Note: As soon as you have selected one of the two dates, we will propose directly the earliest date accordingly to our custom product availability times. Dates shown in pale grey are dates on which delivery of personalized products can only be guaranteed if the date is changed on the opposite calendar. Dates on which delivery is impossible will be indicated in dark grey and disabled.

Cancel

Confirm

Suspend or cancel a case

Display the case details to suspend/cancel it

Click on “suspend surgery” or “cancel surgery”

Waiting for CT scan

DOE John **LEFT KNEE** **[126378]**

[← Back to surgical cases](#)

Medical information

Joint to operate	Knee
Side	Left
Case ID	126378
Preop. extension	5° Flessum
Preop. flexion	30°
Preop. deformity	Varus 12°
Passive correction	Total

[✎ Edit](#)

[⏸ Suspend surgery](#)

[⏹ Cancel surgery](#)

Products to order

Planned implants	ORIGIN CR Femur Cemented ORIGIN CR Fixed Tibia Monobloc Cemented ORIGIN CR Fixed Insert ORIGIN Patella Cemented
Ordered products	ORIGIN CR Femur Cemented ORIGIN CR Fixed Tibia Monobloc Cemented ORIGIN CR Fixed Insert ORIGIN Patella Cemented ORIGIN KNEE-PLAN Femoral Cut Guide ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

[✎ Edit](#)

CT scan exam and surgery date

Radiology center	Wigan Albert Edward Infirmary - WN1 2NN Wigan
CT exam date	09.08.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge

[✎ Edit](#)

Suspend a case (can be resumed later)

Suspend surgery?

By choosing "Yes", you will confirm that you want to suspend the surgical procedure for an undetermined period of time. The surgery date will be reset and the current case will now be classified in 'suspended cases'. Are you sure that you want to suspend this surgical case?

Please enter additional comments if necessary

Cancel Yes

Enter a message / comment to explain reasons why you would like to suspend and click on "yes"

Symbios Connect

Pr CHARNLEY John

Pending planning approval

← Back to surgical cases

⌚ Suspend surgery

⊗ Cancel surgery

2 days remaining to approve planning

📄 Download latest planning.

✅ Approve planning

✖ Modify planning

❓ Ask a question to engineer

🕒 View planning history

EINSTEIN Albert

Patient information

Administrative

Date of birth 18.04.1955 Side Right

Gender Male Case ID 126369

Comments

ORIGIN®

products

- ORIGIN CR Fixed Tibia Monobloc Cemented
- ORIGIN CR Fixed Insert
- ORIGIN Patella Cemented
- ORIGIN KNEE-PLAN Femoral Cut Guide
- ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

Display the suspended case(s) in the list



Symbios Con

Deactivate all status filters and activate the “suspended” filter to easily find the case

Pr CHARNLEY John ▾

Surgical cases

New case

Search



25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
-	CHANEL Coco	HIP	L	INDIVIDUAL HIP™	Suspended

Status quick filters

Active cases [18]

- Pending approvals (4)
- In process (13)
- Delayed (1)

Inactive cases [2]

- Operated (1)
- Cancelled (0)
- Suspended (1)
- To complete (0)

Rows / page: 15 30 50

CHANEL Coco
 HIP LEFT (126368) ✕

Related order

INDIVIDUAL HIP™

Personalized hip femoral stem

Case ID 126368

Planned implants INDIVIDUAL HIP Cemented
HILOCK REV
INLOCK X

Ordered products INDIVIDUAL HIP Cemented
INDIVIDUAL HIP Rasp Impaction hole without thread

Resume a suspended case

 CHANEL Coco LEFT HIP [126368]

Suspended

Patient information

Go in case details and click on “resume surgery”

Comments

[← Back to surgical cases](#)

Resume surgery

Cancel surgery

Weight 65 kg
Size 164 cm

Products to order



INDIVIDUAL HIP™

Planned implants	INDIVIDUAL HIP Cemented HILOCK REV INLOCK X
Ordered products	INDIVIDUAL HIP Cemented INDIVIDUAL HIP Rasp Impaction hole without thread

[↻ View planning history](#)

CT scan exam and surgery date

Radiology center	Imagerie Médicale Hospices Civils de Beaune - 21200 BEAUNE
CT exam date	01.07.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date	None

 [Download CT scan](#)

Resume a suspended case – new surgery date

CT scan exam and surgery date

Select medical information

Imagerie Médicale Hos
BEAUNE
avenue Guigone de Salins
21200 BEAUNE
FRANCE

Tel: 03.80.24.44.32
Fax:

Theatre Stores
WN6 9EP Appley Bridge
UNITED KINGDOM

Tel: +44 1942 244000
Fax:

Select a new date for the surgery

As a consequence of suspending the case, the surgery date has been reset. Please select a new surgery date for this case.

OK

As the case was suspended, the surgery date is not anymore defined. Indicate a new surgery date using the calendar that displays the delivery constraints.

10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

It may not be possible to resume the initial surgery date as time may have elapsed and there is still some constraints of production + delivery

Cancel Confirm

History tracking

Resume a suspended case – new surgery date

Symbios Connect

Pr CHARNLEY John ▾

Confirm changes or make other changes

You have the option of confirming the changes made or making other changes at further stages in the order.
Please note: Changing some information may reset other stages in the order, since the items of information are linked.

Other changes Confirm

CHANEL Coco

Patient information

Administrative information	Medical information
Patient name CHANEL Coco	Joint to operate Hip
Date of birth 18.05.1956	Side Left
Gender Female	Case ID 126368

← [Back to surgical cases](#)

[Resume surgery](#)

[Cancel surgery](#)

[View planning history](#)

Planned implants

- INDIVIDUAL HIP Cemented HILOCK REV INLOCK X

Ordered products

- INDIVIDUAL HIP Cemented
- INDIVIDUAL HIP Rasp Impaction hole without thread

INDIVIDUAL HIP™

CT scan exam and surgery date

Radiology center Imagerie Médicale Hospices Civils de Beaune - 21200 BEAUNE

Click on “confirm” to finalize the resume operation

Cancel a case (definitive cancellation)

Cancel surgery?

By choosing 'Yes', you will confirm that you want to definitively cancel the surgical procedure. The surgery date will be reset and the current case will now be classified under 'cancelled cases'. Are you sure that you want to cancel this surgical case?

Please enter additional comments if necessary

Cancel Yes

Enter a message / comment to explain reasons why you would like to cancel and click on "yes"

Symbios Connect

Pr CHARNLEY John ▾

EINSTEIN Albert

Patient information

Administrative

Patient name	EINSTEIN Albert	Joint to operate	Knee
Date of birth	18.04.1955	Side	Right
Gender	Male	Case ID	126369

Comments

ORIGIN CR Fixed Insert
ORIGIN Patella Cemented
ORIGIN KNEE-PLAN Femoral Cut Guide
ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD

Pending planning approval

← Back to surgical cases

- ⏸ Suspend surgery
- ✖ Cancel surgery

🕒 **2 days** remaining to approve planning

- 📄 Download latest planning
- ✅ Approve planning
- ✖ Modify planning
- ❓ Ask a question to engineer
- 🕒 View planning history

HIP-PLAN[®] functions

HIP-PLAN® cases are to facilitate CT / planning sharing



In Symbios Connect, HIP-PLAN® cases are used mainly to:

1. Register the case to allow CT transfer from Symbios Box
2. Share planning report(s) between sales rep and surgeon

CT scan available

[← Back to surgical cases](#)

- [✎ Edit](#)
- [⏸ Suspend surgery](#)
- [⏹ Cancel surgery](#)

Products to order



HIP-PLAN®

Planned implants	SPS Evolution APRIL Ceramic BIOLOX Delta Insert
Ordered products	None

- [✎ Edit](#)
- [📄 Upload planning report](#)

CT scan exam and surgery date

Radiology center	Centre De Radiologie De La Clinique Montreal Le Tivoli - 11000 CARCASSONNE
CT exam date	01.07.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date	30.08.2022

- [✎ Edit](#)
- [📄 Download CT scan](#)

Download CT scan DICOM file

 ROOSEVELT Franklin LEFT HIP [126362]

CT scan available

Patient information



Click on “Download CT scan” to download

[← Back to surgical cases](#)

[✎ Edit](#)

[⏸ Suspend surgery](#)

[⏹ Cancel surgery](#)

Comments

Weight 78 kg
Size 175 cm

Products to order



Planned implants

SPS Evolution
APRIL Ceramic
BIOLOX Delta Insert

Ordered products

None

HIP-PLAN®

[✎ Edit](#)

[📄 Upload planning report](#)

CT scan exam and surgery date

Radiology center	Centre De Radiologie De La Clinique Montreal Le Tivoli - 11000 CARCASSONNE
CT exam date	01.07.2022
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date	30.08.2022

[✎ Edit](#)

[📄 Download CT scan](#)

Symbios Connect - Version 3.4.3

Share HIP-PLAN® PDF report with surgeon



 ROOSEVELT Franklin LEFT HIP [126362]

CT scan available

Click on “Upload planning report” to upload a PDF HIP-PLAN® report on Symbios server. Other users will be able to download it from Symbios Connect.

[← Back to surgical cases](#)

- [✎ Edit](#)
- [⏸ Suspend surgery](#)
- [⏹ Cancel surgery](#)

Products to order



HIP-PLAN®

Planned implants
SPS Evolution
APRIL Ceramic
BIOLOX Delta Insert

Ordered products
None

- [✎ Edit](#)
- [📄 Upload planning report](#)

CT scan exam and surgery date

Radiology center Centre De Radiologie De La Clinique Montreal Le Tivoli - 11000 CARCASSONNE
CT exam date 01.07.2022

Hospital WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date 30.08.2022

- [✎ Edit](#)
- [📄 Download CT scan](#)

Download HIP-PLAN® PDF report



 ROOSEVELT Franklin LEFT HIP [126362]

Planning available

Click on “Download planning report” to download

Planned implants	SPS Evolution APRIL Ceramic BILOX Delta Insert	
Ordered products	None	



HIP-PLAN®

Radiology center	Centre De Radiologie De La Clinique Montreal Le Tivoli - 11000 CARCASSONNE	
CT exam date	01.07.2022	
Hospital	WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge	
Surgery date	30.08.2022	

[← Back to surgical cases](#)

-  [Edit](#)
-  [Suspend surgery](#)
-  [Cancel surgery](#)

-  [Edit](#)
-  [Upload planning report](#)
-  [Download planning report\(s\)](#)
-  [Ask a question to planner](#)
-  [View planning history](#)

-  [Edit](#)
-  [Download CT scan](#)

Ask question to the planner

 ROOSEVELT Franklin LEFT HIP [126362]

Planning available

Patient info



Click on "Ask a question to planner"

[← Back to surgical cases](#)

Comments

Hip Left
126362

- [✎ Edit](#)
- [⌚ Suspend surgery](#)
- [⊘ Cancel surgery](#)

Weight 78 kg
Size 175 cm

Products to order



HIP-PLAN®

Planned implants SPS Evolution
APRIL Ceramic
BIOLOX Delta Insert

Ordered products None

- [✎ Edit](#)
- [📄 Upload planning report](#)
- [📄 Download planning report\(s\)](#)
- [? Ask a question to planner](#)
- [🕒 View planning history](#)


CT scan exam and surgery date

Radiology center Centre De Radiologie De La Clinique Montreal Le Tivoli - 11000 CARCASSONNE
CT exam date 01.07.2022

Hospital WRIGHTINGTON HOSPITAL - WN6 9EP Appley Bridge
Surgery date 30.08.2022

- [✎ Edit](#)
- [📄 Download CT scan](#)

Ask question to the planner



Ask a question to planner

Please select the planning you would like to ask a particular question about.

Type here any question you want to ask the user who uploaded this planning. The message will be sent to them by email and via Symbios Connect.

Enter your question here

Cancel Submit

To ask a question to planner, proceed as follows:

1. Select which planning you refer to
2. Type question in textbox
3. Click on “submit” button to send question

User preferences

Display user menu - Click on Preferences



Click here to display menus. Click on "preferences".

Pr CHARNLEY John ▾

Pr CHARNLEY John
marketing@symbios.ch

Preferences

Downloads

Help

Terms of Use

Sign out

Surgical cases

New case

« August 2022 »						
Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4
Today						

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]
- Suspended [0]

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN@	Shipped
	HARI Mata	KNEE	R	ORIGIN@	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN@	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN@	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN@	Planning approved
	MONROE Marilyn	KNEE	R	ORIGIN@	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN@	Pending planning approval
	MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available

< Previous 1 2 Next >

Edit account preferences

Symbios Connect Pr CHARNLEY John ▾

Preferences

- My account**
- Notifications
- Hip planning preferences
- Knee planning preferences

My account

[← Back to surgical cases](#)

Personal information

Title

Last name

First name

Email

Telephone

Website

 [Change my SC password](#)

Display settings

Language

Date format



Edit directly, no need to click on save

Edit notification preferences

Symbios Connect

Pr CHARNLEY John ▾

Preferences

- My account
- Notifications**
- Hip planning preferences
- Knee planning preferences

Notification preferences

- Send me an email when a new case that I will operate on is created
 - By myself
 - By my secretary(ies)
 - By any other user of my Symbios Connect network
- Send me an email when one of my cases is modified, suspended or cancelled
 - By myself
 - By my secretary(ies)
 - By any other user of my Symbios Connect network
- Send me an email when one of my cases changes its status
 - CT available
 - For HIP-PLAN cases
 - For INDIVIDUAL HIP™ cases
 - Planning available for HIP-PLAN cases
 - Pending planning approval
 - Planning approved
 - Manufacturing in progress
 - Ready for delivery

[← Back to surgical cases](#)



Usually most surgeons prefer to avoid being overwhelmed by notification emails – adapt accordingly

Edit hip planning preferences – standard implants

Symbios Connect Pr CHARNLEY John ▾

Preferences

- My account
- Notifications
- Hip planning preferences**
- Standard implants**
- INDIVIDUAL HIP™
- Knee planning preferences

Set standard hip planning preferences in order to automatically select them for HIP-PLAN® cases

[← Back to surgical cases](#)

Hip standard implants preferences

Get faster in your regular practice by selecting default implants in this section: it will allow the system to select your preferred implants automatically during the 'New case' sequence for a HIP-PLAN case. However, you will always have the possibility to change your implant selection for a given case if you want to.

Please select the femoral components that you plan the most frequently in your practice:



SPS Evolution
Cementless anatomic femoral stem. Titanium alloy (Ti6Al4V-ISO 5832-3). Coatings: Porous titanium and hydroxyapatite. 12/14 5°40' taper. CCD shaft angle 129°.



BIOLOX Delta Head

- Change
- None
- SPS Evolution**
- SPS HA
- HARMONY
- HARMONY Cemented
- ARCAD HA
- ARCAD LONG
- ARCAD SO
- ARCAD XL
- ATLANTIS Modular

Edit hip planning preferences – INDIVIDUAL HIP[®] stem

Symbios Conn  Pr CHARNLEY John ▾

Set custom stem preferences in order to automatically select them for INDIVIDUAL HIP[®] cases

Preferences

- My account
- Notifications
- Hip planning preferences**
- Standard implants
- INDIVIDUAL HIP[™]**
- Knee planning preferences

INDIVIDUAL HIP[™] preferences

Get faster in your regular practice by selecting default options in this section: it will allow the system to select your preferred implants automatically during the 'New case' sequence for a INDIVIDUAL HIP[™] case. However, you will always have the possibility to change your implant selection for a given case if you want to.

Please select the femoral components that you plan the most frequently in your practice:

[← Back to surgical cases](#)



INDIVIDUAL HIP TiHA

Personalized cementless stem. Titanium alloy (Ti6Al4V-ISO 5832-3), coated with porous titanium and hydroxyapatite. 12/14 5°40' taper. Change ▾

Ref: 3000 0201

Design options

Coating	Collar
<input type="radio"/> Partial HA coating	<input checked="" type="radio"/> No collar
<input checked="" type="radio"/> Full HA coating	<input type="radio"/> Normal collar

Edit hip planning preferences – INDIVIDUAL HIP[®] rasp

Symbios Connect 

Pr CHARNLEY John ▾

Set custom rasp preferences in order to automatically select them for INDIVIDUAL HIP[®] cases

Preferences

- My account
- Notifications
- Hip planning preferences**
 - Standard implants
 - INDIVIDUAL HIP[™]**
 - Knee planning preferences

[← Back to surgical cases](#)

INDIVIDUAL HIP Rasp

Personalized piqued rasp. Stainless steel [316L - NF S 94-090], 12/14 5°40' taper.
Ref: 3000 0701

Design options

Impaction hole

- With thread (use with INDIVIDUAL HIP[™] rasp handle)
- Without thread (use with INDIVIDUAL HIP[™] stem impactor)

[Change](#)

- INDIVIDUAL HIP Compactor
- INDIVIDUAL HIP Compactor NT
- INDIVIDUAL HIP Rasp**
- INDIVIDUAL HIP Rasp NT



BIOLOX Delta Head

Ceramic head [Al2O3 + ZrO2-ISO 6474-2], compatible with 12/14 5°40' taper.

[Change](#)



Edit knee planning preferences – standard implants

Symbios Con 

Pr CHARNLEY John ▾

Set implant preferences in order to automatically select them for FIRST KNEE-PLAN® cases

Preferences

- My account
- Notifications
- Hip planning preferences
- Knee planning preferences**
 - Standard implants**
 - ORIGIN

Knee standard implants preferences

Get faster in your regular practice by selecting default options in this section: it will allow the system to select your preferred implants automatically during the 'New case' sequence for a KNEE-PLAN case. However, you will always have the possibility to change your implant selection for a given case if you want to.

Please select the implants that you plan the most frequently in your practice:

	FIRST Femur Cemented Cemented femoral component. Cobalt-Chrome [CoCrMo-ISO 5832-4].	Change ▾
	FIRST Fixed Tibia Monobloc Cemented Cemented monobloc tibial component for fixed platform. Titanium alloy [Ti6Al4V-ISO 5832-3].	Change ▾

[← Back to surgical cases](#)

Edit knee planning preferences – ORIGIN[®] implants

Symbios C 

Pr CHARNLEY John ▾

Set implant preferences in order to automatically select them for ORIGIN[®] cases

Preferences

- My account
- Notifications
- Hip planning preferences
- Knee planning preferences**
- Standard implants
- ORIGIN**

ORIGIN preferences

Get faster in your regular practice by selecting default options in this section: it will allow the system to select your preferred implants automatically during the 'New case' sequence for an ORIGIN case. However, you will always have the possibility to change your implant selection for a given case if you want to.

Please select the implants that you plan the most frequently in your practice:

	ORIGIN CR Femur Cemented Custom-made cemented femoral component	Change ▾
	ORIGIN CR Fixed Tibia Monobloc Cemented Custom-made cemented monobloc tibial component for fixed platform	Change ▾

[← Back to surgical cases](#)

Edit knee planning preferences – ORIGIN[®] spacer

Symbios Connect Pr CHARNLEY John ▾

Set spacer preference in order to automatically add it to product in future ORIGIN[®] orders

Preferences

- My account
- Notifications
- Hip planning preferences
- Knee planning preferences**
- Standard implants
- ORIGIN**

[← Back to surgical cases](#)



ORIGIN KNEE-PLAN Tibial Cut Guide - INT STD
Custom-made tibial cut guide for internal approaches. Polyamide PA2200.
Delivered sterile with tibial bone model

Please select if you would like the ORIGIN[®] spacer included in the ORIGIN[®] PS Femur Set

 **ORIGIN Spacer**

✓ Change
None
ORIGIN Spacer

Downloads

Display user menu - Click on Downloads



Click here to display menus. Click on "Downloads"

Pr CHARNLEY John ▾

Pr CHARNLEY John
marketing@symbios.ch

- Preferences
- Downloads**
- Help
- Terms of Use
- Sign out

Search

Surgical cases

New case

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]
- Suspended [0]

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN@	Shipped
	HARI Mata	KNEE	R	ORIGIN@	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN@	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN@	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN@	Planning approved
	MONROE Marilyn	KNEE	R	ORIGIN@	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN@	Pending planning approval
	MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available

< Previous 1 2 Next >

Download HIP-PLAN[®] software

Symbios Connect Pr CHARNLEY John ▾

Downloads

	Operating System	Version		System requirements
HIP-PLAN 3	MacOS	3.0.2	Download	MacOSX 10.10 - 10.15 Minimum memory (RAM): 8Gb Minimum hard drive: 20Gb
Symbios Box Desktop	Windows	3.0.2	Download	Windows 7, 8, 8.1 (64bit) Minimum memory (RAM): 8Gb Minimum hard drive: 20Gb Not supported: Integrated Intel graphic card 32bit architectures

 **Select MacOS or Windows version and click on "Download"**

Download Symbios Box Desktop version

Symbios Connect Pr CHARNLEY John ▾

Downloads

	Operating System	Version		System requirements
HIP-PLAN 3	MacOS	3.0.2	Download	MacOSX 10.10 - 10.15 Minimum memory (RAM): 8Gb Minimum hard drive: 20Gb
Symbios Box Desktop	Windows	3.0.2	Download	Windows 7, 8, 8.1 (64bit) Minimum memory (RAM): 8Gb Minimum hard drive: 20Gb Not supported: Integrated Intel graphic card 32bit architectures

 **Select MacOS or Windows version and click on "Download"**

Help menu

Display user menu - Click on Help



Click here to display menus. Click on "Help"

Pr CHARNLEY John ▾

Pr CHARNLEY John
marketing@symbios.ch

- Preferences
- Downloads
- Help
- Terms of Use
- Sign out

Search

Surgical cases

New case

« August 2022 »

Mo	Tu	We	Th	Fr	Sa	Su
25	26	27	28	29	30	31
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31	1	2	3	4

Today

Status quick filters

Active cases [19]

- Pending approvals [4]
- In process [14]
- Delayed [1]

Inactive cases [1]

- Operated [1]
- Cancelled [0]
- Suspended [0]

Surgery date	Patient	Joint ▾	Side ▾	Related order ▾	Status ▾
23.08.2022	HEPBURN Audrey	HIP	L	INDIVIDUAL HIP™	Shipped
	FORD Henri	KNEE	R	ORIGIN@	Shipped
	HARI Mata	KNEE	R	ORIGIN@	Waiting for CT scan
30.08.2022	DESCARTES René	KNEE	R	ORIGIN@	Manufacturing in progress
	BAKER Joséphine	KNEE	L	ORIGIN@	Ready for delivery
	CHOPIN Frédéric	HIP	L	HIP-PLAN	CT scan available
	ROOSEVELT Franklin	HIP	L	HIP-PLAN	CT scan available
06.09.2022	JEFFERSON Thomas	KNEE	L	ORIGIN@	Planning approved
	MONROE Marilyn	KNEE	R	ORIGIN@	Pending planning approval
13.09.2022	EINSTEIN Albert	KNEE	R	ORIGIN@	Pending planning approval
	MANDELA Nelson	HIP	L	HIP-PLAN	CT scan available

< Previous 1 2 Next >

Display features of current version

Pr CHARNLEY John ▾

Symbios Connect

Help

- New features
- Newsletters

Introduction to
Symbios Connect 3



Display newsletters


Symbios Connect

Help

- New features
- Newsletters**
- Launch of ORIGIN® CR
- HIP-PLAN® 3: New version and problem solving
- Total makeover of HIP-PLAN®
- Launch of INDIVIDUAL HIP
- Launch of ORIGIN® PS Modular
- Total makeover of the ORIGIN® instrumentation
- Update of ORIGIN® instruments

ORIGIN® CR Launch

We are pleased to announce the launch of the ORIGIN® CR prosthesis which completes the ORIGIN® patient-matched offer with a posterior cruciate ligament (PCL) preserved version. It will be available in monobloc and modular (modular tibial extension) version.



The ORIGIN® CR will be available for ordering in Symbios Connect **as of 25 January 2022.**

MKT.NEWS0008V1.0.EN

1/6

Symbios Connect - Version 3.4.3

THANK YOU

SYMBIOS®, ORIGIN®, INDIVIDUAL HIP®, KNEE-PLAN® and HIP-PLAN® are registered trademarks belonging to Symbios Orthopédie S.A, Switzerland.

The information contained in this document is intended exclusively for doctors and is by no means presented for diagnostic use or for treating a specific clinical case. It is therefore not a substitute for a doctor's opinion. The products presented in this document are to be used by and under the responsibility of specialized doctors only.

